# Ferritin Testkit

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# Ferritin test kit

**English** For human medical use only!

for quantitative in vitro determination of ferritin on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0170	Ferritin test kit	32 tests
Order number: ST0171	Ferritin test kit	16 tests
Order number: ST1700	Ferritin control kit	2 x 1 ml (decision level)

**Test kit preparation: Allow single test at least 15 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

The serum ferritin concentration declines very early in the development of iron deficiency. A large number of chronic diseases result in increased serum ferritin concentrations. I.e. chronic infections, chronic inflammatory disorders and numerous types of malignancies, especially lymphomas, leukaemias, breast cancer and neuroblastoma. Increase in serum ferritin concentration also occurs in viral hepatitis or following toxic liver injury because of release of ferritin from damaged liver cells. Serum ferritin concentration is also increased with increases of iron stores, as seen in patients with haemosiderosis or haemochromatosis. Besides the use of ferritin as an iron metabolism parameter, ferritin has also gained importance as a tumour marker for therapeutic drug monitoring and follow-up.

## Method

Measurement of antigen/latex-antibody reaction by the kinetic method.

## Measurement Range

Measurement range: 12 - 500 ng/ml  
 Lower detection limit: 12 ng/ml

## Test Kit

ERS cuvette filled with buffer  
 ERS cap filled with antibodies latex reagent

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. The stability gets limited with opening the vacuum package, when stored at 2 - 8 °C, to 3 months beginning from the date of opening. The maximum stability is set by the expiration date stated on the label.

## Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

## Waste Management

Please refer to local legal requirements.

## Specimen Sample Material

Use fresh serum. If the test can not be carried out on the same day, the serum may be stored at 2 - 8 °C for 48 hours. If stored for a longer period, the sample should be frozen. Alternatively, li-hep plasma and EDTA plasma can be used. For capillary samples use capillary li-hep plasma tubes only.

## Reference Range

Men: 20 - 300 ng/ml  
 Women: 15 - 200 ng/ml

This range is given for orientation only. It is recommended that each laboratory establishes its own reference ranges.

## Quality Control

For internal quality control the Ferritin control kit should be used. Order number: ST1700

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**Allow single test at least 15 minutes to warm up to room temperature before use!**

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1.

**CUBESmart**

1.1: A hand holds a test strip labeled 'Ferritin' and inserts it into the 'EURO lysate' compartment of the CUBESmart device.


1.2: A hand holds a vial with a blue cap and dispenses liquid into the test strip.

1.3: The vial is held over the test strip, and the cap is removed, allowing the liquid to flow into the strip.

1.4: The CUBESmart device screen shows a green play button and a test strip icon, indicating the start of the measurement process.

**CUBESmart**

2.

2.1

3.

**CUBESmart**

3.1: A hand uses a pipette to add liquid to a test strip.

3.2: The test strip is held over the device's opening, and the liquid is dispensed.

3.3: A hand inserts the test strip into the device's cartridge slot.

3.4: The CUBESmart device screen shows 'Insert Cartridge' with a green play button and a red 'X' button, and 'Insert cartridge, close door' with an image of the device door being closed.

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### ATTENTION!

Allow single test at least 15 minutes to warm up to room temperature before use!

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### 2. Sample preparation from primary tube

- 2.1 Aspirate 50 µl sample from primary tube

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### 3. Test processing

- 3.1 Dispense 50 µl sample INTO LIQUID in ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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