# HbA1c Testkit

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Correlation with MonoS (20 patient samples):

$$y(\text{Eurolyser HbA1c}) = 0,9438x(\text{Mono S}) - 0,1444; R^2 = 0,9705;$$

Correlation with Menarini HA-8180 (15 patient samples):

$$y(\text{Eurolyser HbA1c}) = 1,0131x(\text{Menarini}) - 0,0857; R^2 = 0,9851;$$

# HbA1c test kit

**English** For human medical use only!

for quantitative in vitro determination of Haemoglobin A1c in human blood on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0110	HbA1c test kit	24 tests
Order number: ST1100	HbA1c control kit	1 x 0.25 ml (decision level)

**Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

Haemoglobin A1c (HbA1c) is a test for the management of patients with diabetes mellitus. Throughout the circulatory life of the red blood cell, HbA1c is formed continuously by adduction of glucose to the N-terminal of the haemoglobin beta chain. This process which is non enzymatic reflects the average exposure of haemoglobin to glucose over an extended period. Studies have shown that HbA1c on diabetic subjects was 2 - 3 fold over levels found in normal individuals.

## Typical Measurement Range

NGSP: 4.0 - 14.0 % (LOT-depending)  
 IFCC: 20 - 130 mmol/mol (LOT-depending)  
 eAG (estimated Average Glucose, LOT-depending)  
 Unit 1: 69 - 359 mg/dl; Unit 2: 3.8 - 19.8 mmol/l

## Limitations

Should not be used for diagnosis of diabetes mellitus. Results may be inconsistent in patients who have the following conditions: Opiate addiction, lead-poisoning, alcoholism, ingestion of large doses of aspirin. Elevated levels of HbF may lead to underestimation of HbA1c. Uremia does not interfere with HbA1c determination. Labile intermediates (Schiff-base) are not detected and therefore, do not interfere with HbA1c determination. Haemoglobin variants HbA2, HbC and HbS do not interfere with this method. Other very rare variants of haemoglobin (e.g. HbE) have not been assessed.

## Principle

Turbidimetric measurement of the HbA1c concentration at 700 nm wavelength.

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2 - 8 °C. The stability gets limited with opening the vacuum package, when stored at 2 - 8 °C, to 3 months beginning from the date of opening.

## Warnings and Precautions

The reagents contain sodium azide (0.95 g/l) as preservative. DO NOT INGEST! Avoid contact with skin and mucous membranes. Take the necessary precautions for the use of laboratory reagents.

## Waste Management

Please refer to local legal requirements.

## Sample Material

Capillary blood from fingertip (fresh) or venous EDTA whole blood (stable for 1 week if maintained at 2 - 8 °C).

## Reference Ranges

IFCC < 42 mmol/mol for non-diabetic (< 53 mmol/mol for glycemic controlled diabetics);  
 NGSP < 6.00 % for non-diabetic (< 7.0 % for glycemic controlled diabetics);

Values are calculated according the master equations:

$$IFCC = (10,93 * NGSP) - 23,50; eAG(mg/dl) = (28,7 * NGSP) - 46,7;$$

## Quality Control

For internal QC use the HbA1c control kit only. Order no.: ST1100

## Precision and Correlation

Imprecision within-run: N = 10; mean = 5.6 % ngsp; CV = 1.2 %

Imprecision day-to-day (each N = 15):

Low (mean = 4.2 % ngsp): CV = 5.2 %  
 Mid (mean = 7.3 % ngsp): CV = 2.5 %  
 High (mean = 9.4 % ngsp): CV = 2.5 %

Correlation with MonoS (20 patient samples):

$$y(\text{Eurolyser HbA1c}) = 0,9438x(\text{Mono S}) - 0,1444; R^2 = 0,9705;$$

Correlation with Menarini HA-8180 (15 patient samples):

$$y(\text{Eurolyser HbA1c}) = 1,0131x(\text{Menarini}) - 0,0857; R^2 = 0,9851;$$

**Literatur/References: 1. Trivelly, LA Ranney H.M and Lai, HAT New Eng J  
 2. Gonen B. and Rubenstein AHDiabetologia 15,1 (1978)**

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### ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature!

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#### 1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place lysis cup, ERS cuvette and ERS cap into test kit rack
- 1.3 Press measurement button, enter required information using the touch screen

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#### 2. Sample preparation

- 2.1a Aspirate 10 µl fresh capillary blood from finger brick using the end-to-end capillary OR ...
- 2.1b Aspirate 10 µl venous whole blood from EDTA primary tube using the pipette
- 2.2a Put filled end-to-end capillary into the lysis cup OR...
- 2.2b Dispense 10 µl sample into lysis cup (prefilled with lysing buffer)
- 2.3 Tightly close and shake lysis cup vigorously 6 to 8 times (= haemolysate) The haemolysate may be used instantly after shaking and is stable for 72 hours at 2 - 8 °C or 8 hours at room temperature.

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#### 3. Sample processing

- 3.1 Pipette 10 µl lysate into ERS cuvette INTO LIQUID
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS Cartridge into laboratory photometer
- 3.4 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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