# PT (INR) Testkit

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# PT (INR) test kit

**English** For human medical use only!  

for quantitative in vitro determination of prothrombin time on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0180	PT (INR) test kit	32 tests
Order number: ST1800	PT (INR) control kit	2 x 1 ml (N & AK)

 **Test kit preparation: Allow single test at least 30 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

PT (INR) test kit is used for the determination of prothrombin complex-activity in whole blood according to Owren. It is ideal for monitoring of oral anticoagulation with vitamin-K inhibitors like Warfarin or Marcumar. PT according to Owren measures the combined activity of the vitamin-K dependent coagulation factors II, VII and X. In analysis, the sample is mixed with thromboplastin, CaCl<sub>2</sub>, fibrinogen and coagulation factor V. The clotting time is prolonged on the therapeutical level (< 0.4 IU/ml) of factor II, V and X in the sample, but the assay is insensitive to heparin in therapeutical levels up to 1.0 IE/ml.

## Assay Principle

Photometrical clot detection.  
 Measurement range: 1 - 8 INR ISI : ~ 1.00 (Lot pending)

## Stability and Storage

**Unopened test kit:**  
 at 2 - 8 °C: Until the expiration date stated on the label.  
 at 20 - 25 °C: 14 days

**R1 reconstituted (100 µl pre-filled in cuvette):**  
 at 2 - 8 °C: 45 days  
 at 20 - 25 °C: 2 days

## Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

## Waste Management

Please refer to local legal requirements.

## Sample Material

20 µl fresh capillary finger blood OR ... 20 µl well mixed citrated whole blood\* (9 + 1 sodium citrate 3.2 %) - not older than 4 hours! \* values may be slightly increased (~ 10 %) due to the dilution factor of the sodium citrate added to the blood collection system.

## Reference Ranges

Typical oral anticoagulant therapeutic range: 2.0 - 3.5 INR

## Quality Control

For internal quality assurance the PT (INR) control kit should be used; Order number: ST1800; Please make sure „Ctrl.Plasma“ is chosen in the sample type menu when performing a measurement of a control.

## EQUAS

For EQUAS procedures we solely can recommend Plasma Control from Technoclone©. Ask your EQUAS institute!

<b>Precision for Plasma Controls</b>	<b>Precision for Citrated Whole Blood</b>
CV (INR 0.93) < 5.75 %	CV (INR 1.31) < 7.05 %
CV (INR 2.27) < 4.30 %	CV (INR 2.41) < 1.28 %

## Correlation for Citrated Whole blood

N=49; R<sup>2</sup>=0.9334; y=0.9273x + 0.174; x=ACL Top300 plasma INR; y=Eurolyser plasma INR

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### 1. Reagent preparation

- 1.1 Open Reagent 1 (R1) and Diluent (D).
- 1.2 Pour Diluent into Reagent 1.
- 1.3 Close Reagent 1.
- 1.4 Gently rotate mixture and allow to reconstitute at room temperature for a minimum of 30 minutes. Meanwhile repeat mixing at least 3 times to achieve good mixture.

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### 2. ERS cuvette preparation (Repeat steps 16x)

- 2.1 Take out cuvette from test kit box and remove cuvette cap.
- 2.2 Aspirate 100 µl from Reagent 1
- 2.3 Dispense 100 µl into ERS cuvette. Close cuvette by use of cuvette cap.
- 2.4 Put back cuvette into test kit package.

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### 3. Post-Processing

- 3.1 Insert date of mixing into the provided label on the test kit box.
- 3.2 Dispose Reagent 1 (remains) and Diluent vial.
- 3.3 Store filled cuvettes upright only (cap to top).

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## ATTENTION!

Testkit is stable for 45 days (stored at 2 - 8 °C) or 2 days (stored at 20 - 25 °C) after reconstitution.

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
**WARNING: Process the following steps within 5 SECONDS!**

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Allow single test at least 30 minutes to warm up to room temperature! Test can be processed on instruments with firmware version 1.33 or higher only.

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- A Insert the plunger into the capillary opening at the wide, black marker. Then, hold capillary in a sloping position to blood drop until it is completely filled up to the small, red marker.  
**Aspirate**
- B WRONG: Capillary is not completely filled with blood up to the small, red marker. Hold the capillary horizontally when aspirating to ensure it gets filled completely.
- C CORRECT: Capillary is completely filled with blood up to the small, red marker.  
**Dispense**
- D WRONG: Plunger has not been pushed all the way into the capillary. Blood is left in capillary.
- E CORRECT: Plunger has been pushed completely into the capillary. White stopper is pressed to the end of the capillary. No blood is left.

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- Place RFID card
- Place ERS cuvette in test kit rack
- Place ERS cap in test kit rack
- Press „Measurement“ button, required information is entered by use of the touch display of the analyser

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### 2. Sample preparation for whole blood

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- Choose sample type: Please make sure „Capil.Blood“ is selected when performing a measurement from capillary finger blood or from citrated whole blood. When performing a measurement from control plasma please make sure „Ctrl. Plasma“ is selected.
- Aspirate 20 µl whole blood from fingertip (See „Correct use of capillary“!)

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- Attention: Entire reagent liquid has to be at the bottom of the cuvette. By pressing the plunger: Dispense 20 µl whole blood into ERS cuvette INTO LIQUID!
- Apply ERS Cap
- Place ERS cartridge into laboratory photometer
- Start automatic sample processing by pressing the  start button on the smart laboratory photometer or by closing the door of the CUBE laboratory photometer.

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