NADAL® Strep A Test (test cassette)

REF 222901N-20

de	Gebrauchsanweisung	2	cs	Návod k použití	37
en	Instructions for use	7	fi	Käyttöohje	42
fr	Instructions d'utilisation	12	sv	Användarinstruktioner	47
es	Instrucciones de uso	17	da	Brugervejledning	52
it	Istruzioni per l'uso	22	no	Bruksanvisning	57
pl	Sposób użycia	27	nl	Gebruiksaanwijzing	61
pt	Instruções de Utilização	32		Symbols	66
				Our Teams	68

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Telefon: 0800 - 10 10 871 seit 1921

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15. Referenzen

1. Levinson ML, Frank PF. Differentiation of group A from other beta hemolytic streptococci with bacitracin. J Bacteriol. 1955 Mar;69(3):284-7.
2. Edwards EA, Phillips IA, Suiter WC. Diagnosis of group A streptococcal infections directly from throat secretions. J Clin Microbiol. 1982 Mar; 15(3): 481-3.
3. Lauer BA, Reller LB, Mirrett S. Effect of atmosphere and duration of incubation on primary isolation of group A streptococci from throat cultures. J Clin Microbiol. 1983 Feb; 17(2): 338-40.
4. Gerber MA, Shulman ST. Rapid diagnosis of pharyngitis caused by group A streptococci. Clin Microbiol Rev. 2004 Jul;17(3):571-80.
5. Wessels MR. Clinical practice. Streptococcal pharyngitis. N Engl J Med. 2011 Feb 17;364(7):648-55.
6. Choby BA. Diagnosis and treatment of streptococcal pharyngitis. Am Fam Physician. 2009 Mar 1;79(5):383-90

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1. Intended Use

The NADAL® Strep A Test is a lateral flow chromatographic immunoassay for the qualitative detection of group A *Streptococcus* antigens in human oropharyngeal swab specimens. The test is intended for use as an aid in the diagnosis of *Streptococcus* A (Strep A) infections in patients showing typical symptoms of a local throat infection (see section 12 'Limitations'). The test procedure is not automated and requires no special training or qualification. The NADAL® Strep A Test is designed for professional use only.

2. Introduction and Clinical Significance

Streptococcus pyogenes belongs to the β -haemolytic group A *Streptococcus* and is a major cause of throat infections such as pharyngitis, tonsillitis and scarlet fever. Symptoms range from a mild sore throat to massive inflammatory reactions, which may be characterised by a severe sore throat, redness and swelling of the mucous membranes, purulent exudate and fever. Scarlet fever is caused by toxin-producing strains. Here, additional symptoms include the typical skin rash, perioral pallor and raspberry tongue. The early diagnosis and treatment of group A streptococcal pharyngitis may reduce the severity of symptoms and the likelihood of later complications, such as rheumatic fever and glomerulonephritis.

Conventional methods for the detection of a group A streptococcal infection are dependent on the isolation and subsequent identification of the pathogen and often require 24-48 hours. Developments in immunological methods for detecting group A *Streptococcus* antigens directly from oropharyngeal swabs support physicians in diagnosing group A streptococcal infections without delay.

Expected values

Pharyngitis can be caused by viruses or bacteria. Pharyngitis caused by group A *Streptococcus* occurs in 5 to 20% of adults, whereas in children it is observed, with corresponding symptoms, in approximately 15 to 30% of cases. Children between the ages of 5 and 15 are the most commonly affected. Although *Streptococcus pyogenes* is common around the world, throat infections are particularly prevalent in temperate climates, occurring most often in the late winter and early spring months.

3. Test Principle

The NADAL® Strep A Test enables the detection of group A *Streptococcus* antigens through the visual interpretation of colour development on the internal test strip. Polyclonal antibodies specific to group A *Streptococcus* antigens are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with the polyclonal antibodies to group A *Streptococcus* antigens which are conjugated to coloured particles and precoated onto the conjugate pad of the test cassette. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there are a sufficient number of group A *Streptococcus* antigens in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, while its absence indicates a negative result. The formation of a coloured line in the control line region (C) serves as a

procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® Strep A test cassettes
- 20 extraction tubes, incl. dropper caps
- Additional material provided according to 93/42/EEC:
Due to possible supply shortages of accessory medical products, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.

a) 20 sterile swabs, CE 0197

Jiangsu Changfeng Medical Industry Co. Ltd.,
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (authorised EU
representative: Llins Service & Consulting GmbH,
Obere Seegasse 34/2, 69124 Heidelberg,
Germany)

b) 20 sterile swabs, CE 0197

CITOTEST LABWARE MANUFACTURING Co. Ltd.,
No.48, Xinxu Road, Haimen, Jiangsu province
(authorised EU representative: WellKang Ltd,
Enterprise Hub, NW Business Complex,
1 Beraghmore Rd., Derry, BT48 8SE, Northern
Ireland)

c) 20 sterile swabs, CE 0197

Hangzhou Yiguoren Biotechnology Co. Ltd., Room
402, Building 2, No. 2628, Yuhangtang Road,
Cangqian Street, Yuhang District, Hangzhou,
311121 Zhejiang, China (authorised EU
representative: Zoustech S.L., Paseo de la
Castellana 141, 28049 Madrid, Spain)

d) 20 sterile swabs, CE 2797

Puritan Medical Products Company LLC,
31 School Street, Guilford, Maine 04443-0149
USA (authorised EU representative: EMERGO
EUROPE, The Hague, The Netherlands)

- 1 bottle of Reagent 1 (white cap)*: 1.0 M sodium nitrite (8 mL), light reddish colour due to pH indicator

Danger
H301: Toxic if swallowed
(see section 7 'Warnings and Precautions')

- 1 bottle of Reagent 2 (red cap)*: 0.4 M acetic acid (8 mL)
- 1 bottle of Positive Control (blue cap)*: non-viable group A *Streptococcus* (0.5 mL)
- 1 reagent holder
- 1 package insert

*containing the following preservative: ProClin™ 300: <0.03%.

No hazard labelling for ProClin™ 300 is required according to Regulation (EC) N° 1272/2008 CLP. Concentrations are below the exemption threshold of <0.03%.

5. Additional Materials Required

- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the entire instructions for use prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not swap caps between bottles.
- Do not use Reagents 1 & 2 if they are discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Do not use swabs from damaged pouches.
- Reagent 1 contains 1.0 M sodium nitrite, which is toxic. If the solution is swallowed, immediately call a poison control centre or physician. Rinse the mouth with water.
- Reagents 1 & 2 are slightly caustic. Avoid contact with eyes or mucous membranes. In the event of accidental contact, wash thoroughly with water.
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection and Preparation

Gently insert a sterile swab into the pharynx and collect secretions by brushing the swab several times against the reddened posterior pharyngeal wall, both tonsils and other inflamed areas. Avoid touching the tongue, teeth and gums.

Note:

It is recommended that swab specimens be processed as soon as possible after collection. If swabs are not to be processed immediately, they should be placed in a sterile, dry, tightly capped tube or bottle and refrigerated. Do not freeze swabs. Swabs can be stored at room temperature (15-30°C) for up to 8 hours or refrigerated (2-8°C) for up to 72 hours. All specimens should be brought to room temperature (15-30°C) prior to testing.

If a bacterial culture is required, lightly roll the swab on a group A selective blood agar plate before using it in the test. The extraction reagents in the test will kill bacteria on the swabs and make it impossible to grow them in bacterial culture after extraction. Please note that if the same swab is used for both inoculation of a bacterial culture plate and testing using a rapid test, this may lead to reduced sensitivity of the test.

9. Test Procedure

Bring tests, specimens, reagents and/or controls to room temperature (15-30°C) prior to testing.

To avoid cross contamination, do not allow the tips of the reagent bottles to come into contact with specimen material.

1. Place a clean extraction tube labeled with the patient or control identification into the designated area of the reagent holder.
2. Holding the bottles vertically, add 6 drops (approximately 240 µL) of reagent 1 (light reddish) to the extraction tube, followed by 4 drops (approximately 160 µL) of reagent 2 (colorless). Then close the bottles.
3. Mix the solution by carefully swirling the extraction tube.

Note: The addition of reagent 2 to reagent 1 changes the colour of the solution from light reddish to light yellowish.

4. Immediately insert the swab into the extraction tube. Swirl the swab, pressing it against the wall of the extraction tube to extract the antigens contained within. Repeat at least 10 times.
5. Let the solution stand for 1 minute at room temperature.
6. Remove the swab, pressing it firmly against the wall of the tube to release as much liquid as possible. Dispose of the swab in accordance with guidelines for the handling of infectious agents.

7. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

8. Place the test cassette on a clean and level surface.

9. Using the dropper cap provided, transfer 3 drops (approximately 100 µL) of the extracted solution from the extraction tube to the specimen well (S) of the test cassette.

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

10. Start the timer.

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

11. Wait for the coloured line(s) to appear. Read the test result after 5 minutes. Do not interpret the result after more than 10 minutes.

11. Quality Control

The internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance. A positive control containing heat-killed group A *Streptococcus* is provided with each test kit.

Operating procedure for external quality control testing

1. Holding the bottles vertically, add 6 drops of reagent 1 and 4 drops of reagent 2 to an extraction tube.
2. Thoroughly mix the positive control by shaking the bottle vigorously. Add 1 drop of the positive control to the tube.
3. Place a clean, sterile swab into the tube and swirl it at least 10 times. Leave the swab in the extraction tube for 1 minute, then release the liquid from the swab tip by rolling it against the wall of the extraction tube and squeezing the extraction tube as the swab is withdrawn. Remove and dispose of the swab.
4. Continue as described in step 7 of 'Test Procedure'.

If the control does not yield a positive result, do not use the tests with patient specimens. Repeat the quality control test or contact your distributor.

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T). The result indicates that group A *Streptococcus* antigens have been detected in the specimen.

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region (C). No line develops in the test line region (T). The result indicates that no group A *Streptococcus* antigens have been detected.

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

12. Limitations

- The NADAL® Strep A Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of group A *Streptococcus* in human oropharyngeal swab specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the concentration of group A *Streptococcus* antigens can be determined with this qualitative test.
- The accuracy of the test depends on the quality of the swab specimen. False negative results may occur due to improper specimen collection or storage. A negative result may also be obtained from patients at the onset of the disease due to low antigen concentration.
- The NADAL® Strep A Test does not differentiate asymptomatic carriers of group A *Streptococcus* from those with a symptomatic infection.
- Respiratory infections, including pharyngitis, can be caused by streptococci of serogroups other than group A as well as other pathogens. A negative test result for group A *Streptococcus* does not exclude infection with other pathogenic microorganisms. Positive results do not preclude co-infections with other pathogens.
- If clinical signs and symptoms are not consistent with test results, additional testing using other clinical methods (e.g. bacterial culture) is recommended.
- The test detects both viable and non-viable group A *Streptococcus* and may yield a positive result in the absence of viable microorganisms. The use of antibiotics or over-the-counter medication may suppress the growth of group A *Streptococcus* in culture despite the presence of microorganisms detectable by rapid antigen tests. This might lead to discrepant results.

- The NADAL® Strep A Test only detects the presence of group A *Streptococcus* antigens in specimens and should not be used as the sole criterion for a diagnosis.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.
- As recommended by many guidelines, patient background and typical symptoms according to 'Centor Score' or 'McIsaac Score' should be taken into account.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

A correlation study between the NADAL® Strep A Test and conventional bacterial culture was performed. Two oropharyngeal swab specimens were collected from patients exhibiting symptoms of pharyngitis. One of the swabs was then used for the inoculation of bacterial culture and the other was used for testing with the NADAL® Strep A Test.

Group A *Streptococcus* was recorded as present or not present.

The results are presented in the following table:

NADAL® Strep A Test	Bacterial culture			
		Positive	Negative	Total
	Positive	99	5	104
	Negative	3	195	198
	Total	102	200	302

Diagnostic sensitivity: 97.06% (91.71% - 98.99%)*

Diagnostic specificity: 97.5% (94.28% - 98.93%)*

Overall agreement: 97.35% (94.86% - 98.65%)*

*95% confidence interval

Analytical performance

Detection limit

8 different group A *Streptococcus* strains were examined at different levels using the NADAL® Strep A Test. The detection limit of the test is approximately 1.0×10^5 colony forming units/specimen for all tested strains. This indicates that the NADAL® Strep A Test detects various group A *Streptococcus* strains with reliable sensitivity.

Measuring range

No adverse effect on T-line formation (prozone effect) was observed for group A *Streptococcus* concentrations up to 1.0×10^7 colony forming units/specimen. Thus, the measuring range of the test is between at least 1.0×10^5 and 1.0×10^7 colony forming units/specimen.

Analytical specificity

Cross-reactivity study

Cross-reactivity studies with microorganisms likely to be found in the respiratory tract were performed using the NADAL® Strep A Test. Strep A negative swabs were spiked with $10 \mu\text{L}$ of the following microorganisms (1.0×10^7 microorganisms/mL) and exhibited negative results.

Microorganism	Microorganism
<i>Bordetella pertussis</i>	<i>Streptococcus</i> , group B
<i>Moraxella catarrhalis</i>	<i>Streptococcus</i> , group F
<i>Candida albicans</i>	<i>Streptococcus</i> , group G

Microorganism	Microorganism
<i>Corynebacterium diphtheriae</i>	<i>Streptococcus</i> , group C
<i>Enterococcus faecalis</i>	<i>Streptococcus mutans</i>
<i>Haemophilus influenzae</i>	<i>Streptococcus pneumoniae</i>
<i>Klebsiella pneumoniae</i>	<i>Streptococcus sanguinis</i>
<i>Neisseria gonorrhoeae</i>	<i>Pseudomonas aeruginosa</i>
<i>Neisseria meningitidis</i>	<i>Serratia marcescens</i>
<i>Neisseria sicca</i>	<i>Staphylococcus aureus</i>
<i>Neisseria subflava</i>	<i>Staphylococcus epidermidis</i>

Interference study

The following substances, normally present in respiratory specimens or artificially introduced into the respiratory tract, were evaluated at the concentrations listed below and showed no interference with the NADAL® Strep A Test.

Substance	Concentration	Substance	Concentration
4-acetamidophenol	10 mg/mL	Zanamivir	1 mg/mL
Menthol throat lozenges	5% w/v	Benzocaine throat spray	5% v/v
Acetylsalicylic acid	20 mg/mL	Blood, type A	2% v/v
Albuterol	0.083 mg/mL	Blood, type B	2% v/v
Amantadine	500 ng/mL	Blood, type AB	2% v/v
Mometasone	500 ng/mL	Blood, type O	2% v/v
Mouthwash Listerine®	5% v/v	Dexamethasone	10 mg/mL
Mouthwash Lion	5% v/v	Dextromethorphan	10 mg/mL
Ascorbic acid chewable tablets	5% w/v	Diphenhydramine HCl	5 mg/mL
Beclomethasone	500 ng/mL	Ibuprofen	10 mg/mL
Nasal spray	5% v/v	Oseltamivir	500 ng/mL
Oxymetazoline	0.05 mg/mL	Phenylephrine	1 mg/mL
Tobramycin	500 ng/mL	Mucin	1 mg/mL

Precision

Repeatability

Repeatability was established by testing 20 replicates of negative, low, medium and high positive specimens using 3 lots of the NADAL® Strep A tests. >99% of the specimens were correctly identified (20/20 correct tests per concentration, 95% confidence interval: 95.5% - 100%). The NADAL® Strep A Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of negative, low, medium and high positive specimens. Testing was performed by 5 operators using 3 independent NADAL® Strep A test lots at 3 different sites on 5 separate days. >99% of the specimens were correctly identified (375/375 correct

tests per concentration, 95% confidence interval: 99% - 100%).
The NADAL® Strep A Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® Strep A Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. Levinson ML, Frank PF. Differentiation of group A from other beta hemolytic streptococci with bacitracin. J Bacteriol. 1955 Mar;69(3):284-7.
2. Edwards EA, Phillips IA, Sulter WC. Diagnosis of group A streptococcal infections directly from throat secretions. J Clin Microbiol. 1982 Mar; 15(3): 481-3.
3. Lauer BA, Reller LB, Mirrett S. Effect of atmosphere and duration of incubation on primary isolation of group A streptococci from throat cultures. J Clin Microbiol. 1983 Feb; 17(2): 338-40.
4. Gerber MA, Shulman ST. Rapid diagnosis of pharyngitis caused by group A streptococci. Clin Microbiol Rev. 2004 Jul;17(3):571-80.
5. Wessels MR. Clinical practice. Streptococcal pharyngitis. N Engl J Med. 2011 Feb 17;364(7):648-55.
6. Choby BA. Diagnosis and treatment of streptococcal pharyngitis. Am Fam Physician. 2009 Mar 1;79(5):383-90

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- Edwards EA, Phillips IA, Suiter WC. Diagnosis of group A streptococcal infections directly from throat secretions. J Clin Microbiol. 1982 Mar; 15(3): 481-3.
- Lauer BA, Reller LB, Mirrett S. Effect of atmosphere and duration of incubation on primary isolation of group A streptococci from throat cultures. J Clin Microbiol. 1983 Feb; 17(2): 338-40.
- Gerber MA, Shulman ST. Rapid diagnosis of pharyngitis caused by group A streptococci. Clin Microbiol Rev. 2004 Jul; 17(3): 571-80.
- Wessels MR. Clinical practice. Streptococcal pharyngitis. N Engl J Med. 2011 Feb 17; 364(7): 648-55.
- Choby BA. Diagnosis and treatment of streptococcal pharyngitis. Am Fam Physician. 2009 Mar 1; 79(5): 383-90.

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15. Bibliographie

- Levinson ML, Frank PF. Differentiation of group A from other beta hemolytic streptococci with bacitracin. J Bacteriol. 1955 Mar; 69(3): 284-7.

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15. Bibliografia

1. Levinson ML, Frank PF. Differentiation of group A from other beta hemolytic streptococci with bacitracin. J Bacteriol. 1955 Mar;69(3):284-7.
2. Edwards EA, Phillips IA, Sulter WC. Diagnosis of group A streptococcal infections directly from throat secretions. J Clin Microbiol. 1982 Mar; 15(3): 481-3.
3. Lauer BA, Reller LB, Mirrett S. Effect of atmosphere and duration of incubation on primary isolation of group A streptococci from throat cultures. J Clin Microbiol. 1983 Feb; 17(2): 338-40.
4. Gerber MA, Shulman ST. Rapid diagnosis of pharyngitis caused by group A streptococci. Clin Microbiol Rev. 2004 Jul;17(3):571-80.
5. Wessels MR. Clinical practice. Streptococcal pharyngitis. N Engl J Med. 2011 Feb 17;364(7):648-55.
6. Choby BA. Diagnosis and treatment of streptococcal pharyngitis. Am Fam Physician. 2009 Mar 1;79(5):383-90.

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15. Referenser

1. Levinson ML, Frank PF. Differentiation of group A from other beta hemolytic streptococci with bacitracin. J Bacteriol. 1955 Mar;69(3):284-7.
2. Edwards EA, Phillips IA, Suiter WC. Diagnosis of group A streptococcal infections directly from throat secretions. J Clin Microbiol. 1982 Mar; 15(3): 481-3.
3. Lauer BA, Reller LB, Mirrett S. Effect of atmosphere and duration of incubation on primary isolation of group A streptococci from throat cultures. J Clin Microbiol. 1983 Feb;17(2): 338-40.
4. Gerber MA, Shulman ST. Rapid diagnosis of pharyngitis caused by group A streptococci. Clin Microbiol Rev. 2004 Jul;17(3):571-80.
5. Wessels MR. Clinical practice. Streptococcal pharyngitis. N Engl J Med. 2011 Feb 17;364(7):648-55.
6. Choby BA. Diagnosis and treatment of streptococcal pharyngitis. Am Fam Physician. 2009 Mar 1;79(5):383-90

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15. Referanser

- Levinson ML, Frank PF. Differentiation of group A from other beta hemolytic streptococci with bacitracin. J Bacteriol. 1955 Mar;69(3):284-7.
- Edwards EA, Phillips IA, Sulter WC. Diagnosis of group A streptococcal infections directly from throat secretions. J Clin Microbiol. 1982 Mar; 15(3): 481-3.
- Lauer BA, Reller LB, Mirrett S. Effect of atmosphere and duration of incubation on primary isolation of group A streptococci from throat cultures. J Clin Microbiol. 1983 Feb; 17(2): 338-40.
- Gerber MA, Shulman ST. Rapid diagnosis of pharyngitis caused by group A streptococci. Clin Microbiol Rev. 2004 Jul;17(3):571-80.
- Wessels MR. Clinical practice. Streptococcal pharyngitis. N Engl J Med. 2011 Feb 17;364(7):648-55.
- Choby BA. Diagnosis and treatment of streptococcal pharyngitis. Am Fam Physician. 2009 Mar 1;79(5):383-90

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